

# INFANT BASSINET BAS 101

CE



# INFANT BASSINET BAS 101

**Single moulded polycarbonate tub**

**Easily removable and washable foam bed and infant tub**

## Over all Dimension

Length	:	765 mm
Width	:	400 mm
Height	:	1055 mm

## Infant bed (Tub) Dimension

Length	:	765 mm
Width	:	400 mm
Height	:	250 mm

<b>Weight</b>	:	20 KG
<b>Maneuverability</b>	:	4" Castors (with brakes)
<b>Mattress</b>	:	High Quality PU Foam
<b>Tray</b>	:	Can be Rotated up to 90°
<b>Coating</b>	:	Epoxy/Powder coated for scratch resistance & Rust protection

## Recommended environmental specifications:

Operating temperature range	:	10°C to 60°C
Transport / Storage temperature	:	-10°C to 60°C
Humidity range	:	5% to 100% RH
Atmospheric pressure	:	500 to 1060 HPA



Call us toll free :

1 800 425 3503



**Phoenix Medical Systems Private Limited**

No. 32/4, Jawaharlal Nehru Salai, Ekkattuthangal, Chennai 600 032, India.  
Tel.: 91-44-2225 1177, 2225 1178 Fax : 91-44-2225 0194

✉ : sales@pmsind.com

🌐 : phoenixmedicalsyste.ms.com