

# Electronic Baby Weighing Scale BWS101



## SPECIFICATIONS

Accuracy.....5 g

Resolution.....1 g

### Weight

Maximum weight.....15 kg

Minimum weight .....10 g

### Display

Display type.....LCD display

### Temperature

Operating Temperature.....-10°C to 50°C

### Power requirements

Input.....230 V AC, 50 Hz

Battery backup.....6 V, 1.2 Ah

Battery working time.....1 hour (when fully charged)

Battery charging time.....2 hours (approximate charged)

### Dimensions of base

Length.....300 mm (12 inches)

Width .....265 mm (10 inches)

Height.....85 mm (3.5 inches)

### Pan

Length.....510 mm (20 inches)

Width .....300 mm (12 inches)

Height.....85 mm (3.5 inches)

Material.....Transparent acrylic

Cat. No: 7.2.3-16/03/31 May' 2016

Call us toll free :  
1 800 425 3503



**Phoenix Medical Systems Private Limited**

No. 32/4, Jawaharlal Nehru Salai, Ekkattuthangal, Chennai 600 032, India.  
Tel.: 91-44-2225 1177, 2225 1178 Fax : 91-44-2225 0194

✉ : sales@pmsind.com

🌐 : phoenixmedicalsyste.ms.com