

# NEONATAL INTENSIVE CARE INCUBATOR INC 100

ROBUST • RELIABLE • USER-FRIENDLY



# STURDY construction STABLE environment

NEONATAL INTENSIVE CARE INCUBATOR INC 100

The **PHOENIX INC 100** has been specially designed to provide the safest and most stable environment for the critical neonate. With a wide range of accessories, it can be tailored to meet any situation arising in the care of the infant.

The excellent nursing facilities, safe and reliable operation, ease of cleaning and sturdy construction of the INC 100 make it a vital part of the equipment in the hospital.

## FEATURES

### Elbow operated ports

Allow easy access and eliminate cross-infection.

### Large front

Provides access for major procedures and easy removal of baby.

### Externally operated Trendelenburg and Fowler positioning

Up to 8° on both sides.

### Easy-to-pull-out baby tray

### Canopy with unique hinges

Easy to clean inside the incubator.



Self-supporting tube and lead apertures (made from autoclavable silicone)

Easy and safe to insert tubes and cords.

Microprocessor-based control panel with unique blower and heater

Easy to read and touch-operated

The **INC100** is backed by the many years of experience **PHOENIX** has in designing and supplying sophisticated neonatal care equipment worldwide.

**PHOENIX** endeavours to meet your every need. If your practice has any special requirement, call **PHOENIX**. We will be happy to provide you a solution.

CDSCO Registration No. : 000380 / 2020

## Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India

Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,

✉ India : sales@pmsind.com, International: exports@pmsind.com

🌐 www.phoenixmedicalsystems.com