

# Electronic Baby Weighing Scale BWS101



## SPECIFICATIONS

Accuracy.....5 g

Resolution.....1 g

### Weight

Maximum weight.....10 kg

Minimum weight .....10 g

### Display

Display type.....LCD display

### Temperature

Operating Temperature.....-10°C to 50°C

### Power requirements

Input.....230 V AC, 50 Hz

Battery backup.....6 V, 1.2 Ah

Battery working time.....1 hour (when fully charged)

Battery charging time.....2 hours (approximate charged)

### Dimensions of base

Length.....300 mm (12 inches)

Width .....265 mm (10 inches)

Height.....85 mm (3.5 inches)

### Pan

Length.....510 mm (20 inches)

Width .....300 mm (12 inches)

Height.....85 mm (3.5 inches)

Material.....Transparent acrylic

## Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India

Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,

✉ India : sales@pmsind.com, International: exports@pmsind.com

🌐 www.phoenixmedicalsyste.ms.com