

NEONATAL INTENSIVE CARE INCUBATOR INC 200

ROBUST • RELIABLE • USER-FRIENDLY



PHOENIX
Phoenix Medical Systems (P) Ltd.

STURDY construction STABLE environment

The **PHOENIX INC 200** has been specially designed to provide the safest and most stable environment for the critical neonate. With a wide range of accessories, it can be tailored to meet any situation arising in the care of the infant.

This incubator provides the caregiver excellent access to the infant. It is user friendly, and its operation is safe and reliable. These features, the ease with which the INC 200 may be cleaned and its sturdy construction make it a vital part of the equipment in the NICU.

Motorized Trendelenburg positioning



The canopy may be opened easily at the top or from the sides. The side walls open downwards using a unique mechanism that requires relatively small clearances at the sides of the incubator. This space-saving feature also means that it is easier to place the baby in the incubator or to remove it from within.

The bed in the incubator slides out conveniently for placing the baby inside or for removal. The bed can also be rotated simultaneously for intubation or other procedures.



NEONATAL INTENSIVE CARE INCUBATOR INC 200

FEATURES

- Clear all round visibility
- Double-wall canopy to reduce heat loss
- Unique Side doors opening from both sides
- Built-in weighing scale
- Temperature (Skin & Air), Servo Humidity and Servo oxygen controls
- 7" Color touch-screen
- Motorized bed tilting and vertical adjustment from both sides
- Trend display for temperature, humidity and SpO₂
- Convenient drawer
- Rotatable bed
- Pulse oximeter





PHOENIX endeavours to meet your every need. If your practice has any special requirement, call **PHOENIX**. We will be happy to provide you a solution.

CE 0123

CDSO Registration No. : 000380 / 2020

Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India

Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,

✉ India : sales@pmsind.com, International: exports@pmsind.com

🌐 www.phoenixmedicalsystems.com