

LABOUR DELIVERY RECOVERY BED LDR 100

EASY • ACCOMMODATING • SUPPORTIVE



Technical Specification

Maximum height of the bed - 1150 mm
Minimum height of the bed - 610 mm
Overall Length 2130 mm

Overall Width with Siderails - 1020 mm
Mattress Size 900mm x 1950 mm x 125 mm
Maximum Head Angle from Horizontal 70°
Maximum back Angle (In Trendelenburg - Like Position) 7°

Length of the mattress

Head section - 750 mm
Seat section - 545 mm
Leg Rest section - 655 mm
Thickness of the mattress - 125 mm

Sliding leg rest for comfortable care

Standard Accessories

SSCatch basin, Mattress, Calf support
Hand grip, IV stand, Individual castors with brake
hand set, Head board

Optional Accessories

Full Leg Support



CDSO Registration No. : 000387 / 2020

Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India

Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,

✉ India : sales@pmsind.com, International: exports@pmsind.com

🌐 www.phoenixmedicalsystems.com