


# INFANT BASSINET BAS 100



## Features & Benefits




### REMOVABLE L TRAY

To store the infant's belonging and medicines and helps cleaning easy




### EASE OF CLEANING

The bassinet has no curved edges, hence cleaning is easy



### TRANSPARENT ACRYLIC DOORS

The doors are foldable, Aiding the end users during medical procedures



### TILTING OF BED

The bed can be tilted 8° both trendlenberg and reverse trendlenberg to aid the infant during digestion and surgery



### Castors with brake

Aids smooth and noiseless mobility of the newborns



CDSO Registration No. : 000383 / 2020

## Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India

Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,

✉ India : sales@pmsind.com, International: exports@pmsind.com

🌐 www.phoenixmedicalsystems.com