

# LABOUR DELIVERY RECOVERY BED LDR 100

EASY • ACCOMMODATING • SUPPORTIVE



## Technical Specification

Maximum height of the bed - 1150 mm  
Minimum height of the bed - 610 mm  
Overall Length 2130 mm

Overall Width with Siderails - 1020 mm  
Mattress Size 900mm x 1950 mm x 125 mm  
Maximum Head Angle from Horizontal 70°  
Maximum back Angle (In Trendelenburg - Like Position) 7°

## Length of the mattress

Head section - 750 mm  
Seat section - 545 mm  
Leg Rest section - 655 mm  
Thickness of the mattress - 125 mm

Sliding leg rest for comfortable care

## Standard Accessories

SSCatch basin, Mattress, Calf support  
Hand grip, IV stand, Individual castors with brake  
hand set, Head board

## Optional Accessories

Full Leg Support



## Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India

Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,

✉ India : sales@pmsind.com, International: exports@pmsind.com

🌐 www.phoenixmedicalsyste.ms.com